

Cd46 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P15529

Item Number: CM0022441



Introduction

High-quality CD46 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Recognizes human, mouse, and rat CD46. Ideal for complement regulation and host-pathogen interaction research.

[Learn More](#)

Parameter	Value
Product Name	CD46 Rabbit pAb
Remarks/Alias	MCP; TLX; AHUS2; MIC10; TRA2.10; CD46
Species	Human
Gene ID (Human)	4179
Gene ID	4179
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-320 of human CD46 (NP_722548.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15529
Protein Weight	44kDa
Shipping	Ice bag