

Cd66A/Ceacam1 Rabbit Pab For Wb Ihc-P Elisa - P13688

Item Number: CM0022723



Introduction

Rabbit polyclonal antibody targeting human CEACAM1 (CD66a), tested in WB, IHC-P, ELISA, with cross-reactivity to mouse and rat. Explore immune checkpoint regulation, cell adhesion, and cancer research.

[Learn More](#)

Parameter	Specification
Product Name	CD66a/CEACAM1 Rabbit pAb
Remarks/Alias	BGP; BGP1; BGPI; CEACAM1
Species	Human
GeneID (Human)	634
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-430 of human CEACAM1 (NP_001703.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13688
Protein Weight	58kDa
Shipping	Ice bag