


Cdk1/2 Rabbit Polymab® Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P06493/P24941

Item Number: CM0006408



Introduction

CDK1/2 Rabbit PolymAb® is a rabbit recombinant monoclonal antibody targeting CDK1/CDK2, validated for WB, IF-P, IHC-P, and ELISA. Reacts with human and mouse. Ideal for cell cycle and cancer research.

[Learn More](#)

Parameter	Value
Product Name	CDK1/2 Rabbit PolymAb®
Remarks/Alias	CDC2; CDC28A; P34CDC2; CDKN2; p33(CDK2)
Species	Human
Gene ID (Human)	983/1017
Gene ID	983/1017
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 199-298 of human CDK1/2 (P24941)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P06493/P24941
Protein Weight	34kDa
Shipping	Ice bag