


Cdk5Rap3 Rabbit Polyclonal Antibody For Wb And Elisa - Q96Jb5

Item Number: CM0018565



Introduction

Rabbit polyclonal antibody against CDK5RAP3 (C53/IC53/LZAP), validated for WB and ELISA in human samples; targets a 57 kDa ufmylation adaptor and cell cycle regulator; ideal for ribosome recycling and cancer signaling research.

[Learn More](#)

Parameter	Value
Product Name	CDK5RAP3 Rabbit pAb
Remarks/Alias	C53; IC53; LZAP; HSF-27; MST016; PP1553; OK/SW-cl.114; CDK5RAP3
Species	Human
Gene ID (Human)	80279
Gene ID	80279
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human CDK5RAP3 (NP_788276.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q96JB5
Protein Weight	57kDa
Shipping	Ice bag