

Cend1 Rabbit Polyclonal Antibody For Wb, Elisa - Q8N111

Item Number: CM0022273

Introduction

Rabbit polyclonal antibody against human CEND1 (BM88), validated for Western blot and ELISA. Cross-reacts with mouse and rat. Ideal for studying neuronal differentiation and cell cycle exit.

[Learn More](#)

Parameter	Value
Product Name	CEND1 Rabbit pAb
Remarks/Alias	BM88; CEND1
Species	Human
GeneID (Human)	51286
GeneID	51286
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-123 of human CEND1 (NP_057648.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N111
Protein Weight	15kDa
Shipping	Ice bag