


Clusterin Alpha Chain Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P10909

Item Number: CM0017915



Introduction

Rabbit polyclonal antibody targeting human Clusterin alpha chain (CLU/APOJ/CLI). Validated for Western blot, IHC-P, and ELISA; cross-reacts with mouse and rat. Ideal for apoptosis, chaperone, and complement research.

[Learn More](#)

Parameter	Value
Product Name	Clusterin alpha chain Rabbit pAb
Remarks/Alias	CLI; AAG4; APOJ; CLU1; CLU2; KUB1; SGP2; APO-J; SGP-2; SP-40; TRPM2; TRPM-2; NA1/NA2; Clusterin alpha chain
Species	Human
GeneID (Human)	1191
GeneID	1191
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 314-449 of human Clusterin alpha chain (NP_001822.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10909
Protein Weight	52kDa
Shipping	Ice bag