


Cmpk1 Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P30085

Item Number: CM0025776



Introduction

High-quality CMPK1 Rabbit Polyclonal Antibody validated for WB, IF/ICC, IP, and ELISA. Recognizes human and mouse UMP-CMP kinase, involved in pyrimidine biosynthesis. Ideal for nucleotide metabolism research.

[Learn More](#)

Parameter	Value
Product Name	CMPK1 Rabbit pAb
Remarks/Alias	CK; CMK; UMK; CMPK; UMPK; UMP-CMPK; CMPK1
Species	Human
GeneID (Human)	51727
GeneID	51727
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 69-228 of human CMPK1 (NP_057392.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P30085
Protein Weight	22kDa
Shipping	Ice bag