

Ctcf Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa, Chip - P49711

Item Number: CM0016353



Introduction

CTCF Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, IP, ELISA, and ChIP applications. Detects endogenous CTCF in human, mouse, and rat. Ideal for chromatin biology, gene regulation, and 3D genome organization research.

[Learn More](#)

Parameter	Value
Product Name	CTCF Rabbit pAb
Remarks/Alias	MRD21; FAP108; CFAP108; CTCF
Species	Human
GeneID (Human)	10664
GeneID	10664
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human CTCF (NP_006556.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat
SWISS	P49711
Protein Weight	83 kDa
Shipping	Ice bag