


## Dgcr8 Rabbit Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q8WYq5

Item Number: CM0019748



### Introduction

High-quality DGCR8 rabbit polyclonal antibody for WB, IHC-P, IP, ELISA. Recognizes human, mouse, rat DGCR8 (86 kDa), a key microprocessor subunit critical for miRNA biogenesis and DNA repair.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | DGCR8 Rabbit pAb                                                                                                   |
| Remarks/Alias    | Gy1; pasha; DGCRK6; C22orf12; DGCR8                                                                                |
| Species          | Human                                                                                                              |
| GeneID (Human)   | 54487                                                                                                              |
| GeneID           | 54487                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 180-340 of human DGCR8 (NP_073557.3) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, IHC-P, IP, ELISA                                                                                               |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | Q8WYQ5                                                                                                             |
| Protein Weight   | 86kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |