


Dimethyl-Histone H3 (K79) Rabbit Polyclonal Antibody (Q16695/P68431) For Wb, Elisa, Chip

Item Number: CM0000253



Introduction

DiMethyl-Histone H3-K79 Rabbit pAb validated for Western blot, ELISA, and ChIP. Recognizes human, mouse, and rat H3K79me2. Suitable for epigenetic research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	DiMethyl-Histone H3-K79 Rabbit pAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; DiMethyl-Histone H3-K79
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	Synthetic peptide (A synthetic dimethylated peptide around K79 of human histone H3 (NP_003520.1))
Source	Rabbit
Category	Methylated Antibodies
Application	WB, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag