

Dimethyl-Histone H3-K9 Rabbit Monoclonal Antibody (Q16695/P68431) For Wb, Ihc-P, Elisa, Chip

Item Number: CM0000397



Introduction

Affinity-purified rabbit monoclonal antibody against histone H3 dimethylated at lysine 9 (H3K9me2). Suitable for WB, DB, IHC-P, ELISA, and ChIP. Cross-reacts with human, mouse, rat, and a wide range of species. Ideal for epigenetic and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	DiMethyl-Histone H3-K9 Rabbit mAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; DiMethyl-Histone H3-K9
Species	Human
GeneID (human)	8290/8350
GeneID	8290/8350
Immunogen	Synthetic peptide A synthetic dimethylated peptide around K9 of human Histone H3 (NP_003520.1).
Source	Rabbit
Category	Methylated Antibodies
Application	WB, DB, IHC-P, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag