

Dlgap5 Rabbit Pab For Wb, If/Icc, Elisa - Q15398

Item Number: CM0020862



Introduction

High-quality DLGAP5 rabbit polyclonal antibody validated for Western blot, IF/ICC, and ELISA applications. Targets human DLGAP5 (DLG7/HURP), a cell cycle regulator and potential cancer biomarker.

[Learn More](#)

Parameter	Value
Product Name	DLGAP5 Rabbit pAb
Remarks/Alias	DLG7; HURP; DLGAP5
Species	Human
GeneID (Human)	9787
GeneID	9787
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 547-846 of human DLGAP5 (NP_055565.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q15398
Protein Weight	95kDa
Shipping	Ice bag