

Dmp1 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q13316

Item Number: CM0005574

Introduction

Anti-DMP1 rabbit monoclonal antibody validated for WB, IF/ICC, and ELISA. Targets human dentin matrix acidic phosphoprotein 1 (DMP1), 56kDa, involved in bone and tooth mineralization.

[Learn More](#)

Parameter	Value
Product Name	DMP1 Rabbit mAb
Remarks/Alias	ARHP; ARHR; DMP-1
Species	Human
GeneID	1758
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 413-513 of human DMP1 (NP_004398.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q13316
Protein Weight	56kDa
Shipping	Ice bag