


E2F1 Rabbit Polyclonal Antibody For Wb And Elisa - Q01094

Item Number: CM0023645



Introduction

High-quality rabbit polyclonal antibody against human transcription factor E2F1, validated for WB and ELISA. Suitable for human, mouse, and rat samples. Ideal for cell cycle and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	E2F1 Rabbit pAb
Remarks/Alias	RBP3; E2F-1; RBAP1; RBBP3; E2F1
Species	Human
GeneID (human)	1869
GeneID	1869
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 301-437 of human E2F1 (NP_005216.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01094
Protein Weight	47kDa
Shipping	Ice bag