


Eif2S3 Rabbit Polyclonal Antibody For Wb, Elisa - P41091

Item Number: CM0027734



Introduction

Reliable EIF2S3 rabbit polyclonal antibody for WB and ELISA, cross-reactive with rat. Targets translation initiation factor eIF2-gamma, suitable for protein synthesis research and eIF2 complex studies.

[Learn More](#)

Parameter	Value
Product Name	EIF2S3 Rabbit pAb
Remarks/Alias	EIF2; EIF2G; MEHMO; MRXSBRK; eIF-2gA; EIF2gamma; EIF2S3
Species	Human
GeneID (Human)	1968
GeneID	1968
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 173-472 of human EIF2S3 (NP_001406.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P41091
Protein Weight	51kDa
Shipping	Ice bag