

## Eif4A1 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P60842

Item Number: CM0006734



### Introduction

High-quality rabbit monoclonal antibody against human eIF4A1 (DDX2A), validated for WB, IHC-P, IF/ICC, IP, and ELISA. Recognizes eIF4A1 from human, mouse, rat, and monkey—ideal for translation initiation research.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | eIF4A1 Rabbit mAb                                                                                                 |
| Remarks/Alias    | DDX2A; EIF4A; EIF-4A; eIF4A-I; eIF-4A-I; eIF4A1                                                                   |
| Species          | Human                                                                                                             |
| Gene ID (Human)  | 1973                                                                                                              |
| Gene ID          | 1973                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human eIF4A1 (NP_001407.1) |
| Source           | Rabbit                                                                                                            |
| Category         | Monoclonal Antibodies                                                                                             |
| Application      | WB, IHC-P, IF/ICC, IP, ELISA                                                                                      |
| Cross Reactivity | Human, Mouse, Rat, Monkey                                                                                         |
| SWISS            | P60842                                                                                                            |
| Protein Weight   | 46kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |