

Fshb Rabbit Monoclonal Antibody For If-P, Ihc-P, Elisa - P01225

Item Number: CM0010310



Introduction

Validate FSHB rabbit monoclonal antibody for IF-P, IHC-P, ELISA. Cross-reacts with human, mouse, rat FSHB. Ideal for reproductive biology research on follicle development and spermatogenesis, targeting follicle-stimulating hormone beta subunit.

[Learn More](#)

Parameter	Value
Product Name	FSHB Rabbit mAb
Remarks/Alias	HH24
Species	Human
GeneID (Human)	2488
GeneID	2488
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 25-116 of human FSH alpha/beta(NP_000726.1/NP_000501.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01225
Protein Weight	15 kDa
Shipping	Ice bag