

GLI1 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P08151

Item Number: CM0007617



Introduction

High-quality GLI1 rabbit monoclonal antibody validated for Western blot, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat GLI1; immunogen covers N-terminal amino acids 1-233. Ideal for cancer and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	GLI1 Rabbit mAb
Remarks/Alias	GLI; PPD1; PAPA8; GLI1
Species	Human
GeneID (human)	2735
GeneID	2735
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-233 of human GLI1 (NP_005260.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08151
Protein Weight	118kDa
Shipping	Ice bag