


Gmn Rabbit Pab For Wb, Ihc-P, If/Icc, Elisa - 075496

Item Number: CM0028911



Introduction

Rabbit polyclonal antibody targeting human Geminin (GMNN) for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Validated for research use in cell cycle and DNA replication studies.

[Learn More](#)

Parameter	Value
Product Name	GMNN Rabbit pAb
Remarks/Alias	Gem; MGORS6; GMNN
Species	Human
GeneID (Human)	51053
GeneID	51053
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075496
Protein Weight	24kDa
Shipping	Ice bag