

Gpr182 Rabbit Pab For Wb, Elisa - O15218

Item Number: CM0024148

Introduction

Rabbit polyclonal antibody against GPR182 (ACKR5), an atypical chemokine scavenger receptor. Validated for WB and ELISA, cross-reacts with human and mouse. Ideal for chemokine signaling and immunology research.

[Learn More](#)

Parameter	Value
Product Name	GPR182 Rabbit pAb
Remarks/Alias	AMR, 7TMR, ADMR, AM-R, G10D, L1-R, gamrh, hrhAMR, GPR182
Species	Human
GeneID	11318
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 250-350 of human GPR182 (NP_009195.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O15218
Protein Weight	45kDa
Shipping	Ice bag