


Gsta2 Rabbit Polyclonal Antibody For Wb, Elisa - P09210

Item Number: CM0025213



Introduction

Rabbit anti-GSTA2 polyclonal antibody validated for WB and ELISA. Recognizes human, mouse, rat GSTA2. Recombinant protein immunogen covering full-length human GSTA2. Ideal for detoxification and liver studies. - P09210

[Learn More](#)

Parameter	Value
Product Name	GSTA2 Rabbit pAb
Remarks/Alias	GST2; GTA2; GTH2; GSTA2-2; GSTA2
Species	Human
GeneID	2939
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-222 of human GSTA2 (NP_000837.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09210
Protein Weight	26kDa
Shipping	Ice bag