

Gtf2l Rabbit Polyclonal Antibody For Wb, IHC-P, Elisa - P78347

Item Number: CM0029261



Introduction

Rabbit polyclonal antibody against GTF2I for WB, IHC-P, ELISA. Reacts with Human, Mouse, Rat. General transcription factor II-I involved in transcriptional regulation and implicated in Williams-Beuren syndrome.

[Learn More](#)

Parameter	Value
Product Name	GTF2I Rabbit pAb
Remarks/Alias	WBS; DIWS; SPIN; IB291; BAP135; BTKAP1; TFII-I; WBSCR6; GTFII-I; GTF2I
Species	Human
GeneID (Human)	2969
GeneID	2969
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 600-810 of human GTF2I (NP_001157108.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P78347
Protein Weight	112kDa
Shipping	Ice bag