

Hars2 Rabbit Pab For Wb, Elisa - P49590

Item Number: CM0024897

Introduction

HARS2 rabbit polyclonal antibody validated for WB and ELISA. Detects human, mouse, rat HARS2, a mitochondrial histidine--tRNA ligase. Ideal for studying mitochondrial translation and related diseases.

[Learn More](#)

Parameter	Value
Product Name	HARS2 Rabbit pAb
Remarks/Alias	HO3; HARSL; HARSR; HisRS; PRLTS2; HARS2
Species	Human
GeneID	23438
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 217-506 of human HARS2 (NP_036340.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49590
Protein Weight	57kDa
Shipping	Ice bag