


Human Complement C6 Rabbit Polyclonal Antibody For Wb And Elisa - P13671

Item Number: CM0014279



Introduction

High-quality rabbit polyclonal antibody targeting human complement C6, a key MAC component. Validated for WB and ELISA, with cross-reactivity to mouse and rat. Suitable for complement pathway research.

[Learn More](#)

Parameter	Value
Product Name	C6 Rabbit pAb
Remarks/Alias	C6
Species	Human
Gene ID	729
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 586-934 of human C6 (NP_000056.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P13671
Protein Weight	105kDa
Shipping	Ice bag