

Hus1B Rabbit Pab For Wb, Elisa - Q8Nhy5

Item Number: CM0030632

Introduction

HUS1B rabbit polyclonal antibody for WB and ELISA applications. Recognizes human HUS1B protein (~31 kDa) with cross-reactivity to mouse and rat. Suitable for studies in checkpoint biology and tissue expression analysis.

[Learn More](#)

Parameter	Specification
Product Name	HUS1B Rabbit pAb
Remarks/Alias	HUS1B
Species	Human
GeneID (Human)	135458
GeneID	135458
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-278 of human HUS1B (NP_683762.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NHY5
Protein Weight	31kDa
Shipping	Ice bag