

Ikkα Rabbit Monoclonal Antibody For Wb, Elisa - 015111

Item Number: CM0002604



Introduction

IKKα rabbit monoclonal antibody validated for Western blot and ELISA in human and mouse samples. Immunogen corresponds to N-terminal peptide; expected molecular weight 85kDa. Ideal for NF-kappa-B signaling research.

[Learn More](#)

Parameter	Value
Product Name	IKKα Rabbit mAb
Remarks/Alias	BPS2; IKK1; IKKA; IKBKA; TCF16; NFKBIKA; IKK-alpha; IKKα
Species	Human
GeneID (Human)	1147
GeneID	1147
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human IKKα (NP_001269.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	015111
Protein Weight	85kDa
Shipping	Ice bag