


IL-1 Alpha Rabbit Pab For Wb, Ihc-P, Elisa - P01583

Item Number: CM0015903



Introduction

Rabbit polyclonal antibody targeting human IL-1 alpha (31kDa), validated in WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for studying inflammation, innate immunity, and alarmin signaling.

[Learn More](#)

Parameter	Description
Product Name	IL-1 alpha Rabbit pAb
Remarks/Alias	IL1; IL-1A; IL1F1; IL1-ALPHA; IL-1 alpha
Species	Human
GeneID (human)	3552
GeneID	3552
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 113-267 of human IL-1 alpha (NP_000566.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01583
Protein Weight	31kDa
Shipping	Ice bag