


IL18 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q14116

Item Number: CM0006418



Introduction

High-quality Rabbit monoclonal antibody against human IL18 (Interleukin-18). Validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Ideal for Th1/NK cell immunity and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	IL18 Rabbit mAb
Remarks/Alias	IGIF; IL-18; IL-1g; IL1F4; IL18
Species	Human
GeneID (Human)	3606
GeneID	3606
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-193 of human IL18 (NP_001553.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14116
Protein Weight	22kDa
Shipping	Ice bag