

Inppl1 Rabbit Polyclonal Antibody For Wb And Elisa - O15357

Item Number: CM0013309



Introduction

INPPL1 (SHIP2) rabbit polyclonal antibody validated for WB and ELISA; reacts with human, mouse, and rat. Targets the phosphatidylinositol phosphatase SHIP2, a key regulator of PI3K signaling, for metabolic and cancer research.

[Learn More](#)

Parameter	Value
Product Name	INPPL1 Rabbit pAb
Remarks/Alias	OP5MD; SHIP2; INPPL1
Species	Human
Gene ID (Human)	3636
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-310 of human INPPL1 (NP_001558.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15357
Protein Weight	139kDa
Shipping	Ice bag