


Kcnc1 Rabbit Polyclonal Antibody For Wb, Elisa - P48547

Item Number: CM0025424



Introduction

KCNC1 rabbit polyclonal antibody validated for WB and ELISA, with broad reactivity across human, mouse, and rat. Detects voltage-gated potassium channel KCNC1 (Kv3.1), critical for rapid repolarization in fast-firing brain neurons. Suitable for neuroscience and ion channel research.

[Learn More](#)

Parameter	Value
Product Name	KCNC1 Rabbit pAb
Remarks/Alias	KV4; EPM7; NGK2; KV3.1; KCNC1
Species	Human
GeneID (Human)	3746
GeneID	3746
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human KCNC1 (NP_004967.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48547
Protein Weight	58kDa
Shipping	Ice bag