

Kd Validated Anti-Abca6 Rabbit Polyclonal Antibody For Wb And Elisa - Q8N139

Item Number: CM0020664

Introduction

KD-validated rabbit polyclonal antibody targeting human ABCA6. Suitable for WB and ELISA, with reactivity across human, mouse, and rat. ABCA6 is a probable macrophage lipid transporter localized to the Golgi.

[Learn More](#)

Name	Value
Product Name	[KD Validated] ABCA6 Rabbit pAb
Remarks/Alias	EST155051; [KD Validated] ABCA6
Species	Mouse
GeneID (Human)	23460
GeneID	23460
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1518-1617 of human ABCA6 (NP_525023.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N139
Protein Weight	184kDa
Shipping	Ice bag