

Kd Validated Anti-Dnmt1 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - P26358

Item Number: CM0003731



Introduction

KD-validated recombinant rabbit monoclonal antibody targeting human DNMT1. Validated for WB, IP, ELISA. Cross-reacts with human, mouse, rat. Suitable for studying DNA methylation, epigenetic inheritance, and cancer development.

[Learn More](#)

Parameter	Specification
Product Name	[KD Validated] DNMT1 Rabbit mAb
Remarks/Alias	AIM; DNMT; MCMT; CXXC9; HSN1E; ADCADN; m.Hsa1; [KD Validated] DNMT1
Species	Human
Gene ID (Human)	1786
Gene ID	1786
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P26358
Protein Weight	183kDa
Shipping	Ice bag