

[Kd Validated] Anti-Fgfr3 Polyclonal Antibody For Wb, Fc, Elisa - P22607

Item Number: CM0019096



Introduction

KD validated rabbit polyclonal antibody targeting human FGFR3 (P22607), suitable for Western blot, flow cytometry, and ELISA. FGFR3, a receptor tyrosine kinase, plays key roles in skeletal development and oncogenesis, making it a critical research target.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated]FGFR3 Rabbit pAb
Remarks/Alias	FGFR3; ACH; CD333; CEK2; HSFGR3EX; JTK4; fibroblast growth factor receptor 3
Species	Human
GeneID (Human)	2261
GeneID	2261
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-375 of human FGFR3 (NP_000133.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Human
SWISS	P22607
Protein Weight	75kDa/85kDa/87kDa/88kDa
Shipping	Ice bag