

Kd Validated Anti-Hnrnpa2B1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P22626

Item Number: CM0026205



Introduction

KD validated rabbit polyclonal antibody against human HNRNPA2B1. Suitable for WB, IHC-P, ELISA; reacts with human, mouse, rat. Ideal for RNA processing, transport, and innate immunity studies.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] HNRNPA2B1 Rabbit pAb
Remarks/Alias	RNPA2; HNRPA2; HNRPB1; SNRPB1; HNRNPA2; HNRNPB1; IBMPFD2; HNRPA2B1; B1
Species	Human
Gene ID (Human)	3181
Gene ID	3181
Immunogen	Synthetic peptide corresponding to amino acids 150-250 of human HNRNPA2B1 (NP_112533.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22626
Protein Weight	37kDa
Shipping	Ice bag