

[Kd Validated] Anti-Ptbp1 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P26599

Item Number: CM0007173



Introduction

KD Validated anti-PTBP1 rabbit monoclonal antibody validated for WB, IHC-P, IF/ICC, IP, and ELISA applications. Detects human, mouse, and rat PTBP1. Ideal for pre-mRNA splicing and alternative splicing regulation research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] PTBP1 Rabbit mAb
Remarks/Alias	PTB; PTB2; PTB3; PTB4; pPTB; HNRPI; PTB-1; PTB-T; HNRNPI; HNRNP-I; P1
Species	Human
GeneID	5725
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PTBP1 (NP_114368.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P26599
Protein Weight	60kDa
Shipping	Ice bag