

[Kd Validated] Anti-Rarg Polyclonal Antibody For Wb And Elisa - P13631

Item Number: CM0021381



Introduction

KD validated rabbit polyclonal antibody targeting human retinoic acid receptor gamma (RARG). Suitable for Western blot and ELISA. Detects a ~50 kDa protein involved in gene regulation, development, and skeletal growth.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] RARG Rabbit pAb
Remarks/Alias	RARC; NR1B3; RARGgamma; [KD Validated] RARG
Species	Human
Gene ID (human)	5916
Gene ID	5916
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-70 of human RARG (NP_000957.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P13631
Protein Weight	50kDa
Shipping	Ice bag