

[Kd Validated] Anti-Thymidylate Synthase Monoclonal Antibody For Wb, Ihc-P, Elisa - P04818

Item Number: CM0006205



Introduction

Validated recombinant rabbit monoclonal antibody against thymidylate synthase (TYMS), a critical enzyme in nucleotide biosynthesis. Suitable for WB, IHC-P, ELISA; human reactivity. For cancer research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] Thymidylate synthase Rabbit mAb
Remarks/Alias	TS; TMS; DKCD; HST422
Species	Human
Gene ID	7298
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-313 of human Thymidylate synthase (NP_001062.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P04818
Protein Weight	36kDa
Shipping	Ice bag