


[Kd Validated] Mct1 Rabbit Monoclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - P53985

Item Number: CM0004057



Introduction

Rabbit monoclonal antibody targeting MCT1 (SLC16A1) for WB, IF, IHC-P, ELISA. KD validated, detects human monocarboxylate transporter 1, 54kDa. Ideal for metabolism and cancer research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] MCT1 Rabbit mAb
Remarks/Alias	MCT; HHF7; MCT1; MCT1D
Species	Human
Gene ID (Human)	6566
Gene ID	6566
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 190-245 of human MCT1 (NP_003042.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P53985
Protein Weight	54kDa
Shipping	Ice bag