


[Kd Validated] Pin1 Rabbit Mab For Wb And Elisa - Q13526

Item Number: CM0006928



Introduction

KD validated rabbit monoclonal antibody targeting human Pin1 (UniProt Q13526) for Western blot (WB) and ELISA. Cross-reacts with mouse and rat. Detects endogenous Pin1 at 18 kDa.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] Pin1 Rabbit mAb
Remarks/Alias	DOD; UBL5; [KD Validated] Pin1
Species	Human
Gene ID (Human)	5300
Gene ID	5300
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-163 of human Pin1 (NP_006212.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13526
Protein Weight	18kDa
Shipping	Ice bag