

KLK7 Rabbit Polyclonal Antibody For Wb And Elisa - P49862

Item Number: CM0026786



Introduction

Rabbit polyclonal antibody targeting human KLK7 (kallikrein-7), suited for WB and ELISA. Reacts with human and mouse. Ideal for skin desquamation, ovarian cancer, and protease research.

[Learn More](#)

Parameter	Value
Product Name	KLK7 Rabbit pAb
Remarks/Alias	hK7; SCCE; PRSS6; KLK7
Species	Human
Gene ID (Human)	5650
Gene ID	5650
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-253 of human KLK7 (NP_644806.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P49862
Protein Weight	28kDa
Shipping	Ice bag