

[KO Validated] Anti-Atg5 Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q9H1Y0

Item Number: CM0014691



Introduction

Rabbit polyclonal antibody targeting human ATG5 (Q9H1Y0), validated for WB, IF/ICC, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for autophagy, apoptosis, and viral infection research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] ATG5 Rabbit pAb
Remarks/Alias	ASP; APG5; APG5L; hAPG5; SCAR25; APG5-LIKE; ATG5
Species	Human
Gene ID (Human)	9474
Gene ID	9474
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-275 of human ATG5 (NP_004840.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H1Y0
Protein Weight	32 kDa
Shipping	Ice bag