


KO Validated Anti-Carbonic Anhydrase 2 (CA2) Rabbit Polyclonal Antibody For Wb, Elisa - P00918

Item Number: CM0023101



Introduction

KO validated rabbit polyclonal antibody targeting human Carbonic Anhydrase 2 (CA2). Suitable for western blot (WB) and ELISA applications; cross-reacts with human and mouse. For research use in pH regulation, bone biology, and ion transport studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Carbonic Anhydrase 2 (CA2) Rabbit pAb
Remarks/Alias	CAC; CAII; Car2; CA-II; HEL-76; HEL-S-282
Species	Human
Gene ID (Human)	760
Gene ID	760
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human Carbonic Anhydrase 2 (CA2) (NP_000058.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P00918
Protein Weight	29kDa
Shipping	Ice bag