

Ko Validated Anti-Cypor (Por) Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P16435

Item Number: CM0027495



Introduction

KO validated rabbit polyclonal antibody against human NADPH--cytochrome P450 reductase (CYPOR/POR). Suitable for WB, IHC-P, IF/ICC, and ELISA. Recognizes human, mouse, and rat CYPOR. Designed for IVD research and assay development.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CYPOR Rabbit pAb
Remarks/Alias	CPR; CYPOR; P450R; OR
Species	Human
Gene ID (Human)	5447
Gene ID	5447
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human CYPOR (NP_000932.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16435
Protein Weight	77kDa
Shipping	Ice bag