


[KO Validated] Anti-Esd Rabbit Polyclonal Antibody For Wb And Elisa - P10768

Item Number: CM0023032



Introduction

KO validated rabbit polyclonal antibody targeting human ESD (S-formylglutathione hydrolase, esterase D). Cross-reacts with mouse, supports WB and ELISA for formaldehyde metabolism research, ideal for studying detoxification pathways.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] ESD Rabbit pAb
Remarks/Alias	FGH; SD
Species	Human
Gene ID	2098
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-282 of human ESD (NP_001975.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P10768
Protein Weight	31kDa
Shipping	Ice bag