


[KO Validated] Anti-Fam38A/Piezo1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q92508

Item Number: CM0006655



Introduction

KO-validated rabbit monoclonal antibody targeting PIEZO1, a mechanosensitive ion channel involved in touch, vascular development, and cell adhesion. Suitable for WB, IHC-P, and ELISA; detects human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] FAM38A/PIEZO1 Rabbit mAb
Remarks/Alias	ER; DHS; Mib; LMPH3; FAM38A; LMPHM6; FAM38A/PIEZO1
Species	Human
Gene ID (Human)	9780
Gene ID	9780
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2199-2429 of human PIEZO1 (NP_001136336.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92508
Protein Weight	287kDa
Shipping	Ice bag