

KO Validated Anti-Insulin-Degrading Enzyme (Ide) Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P14735

Item Number: CM0022841



Introduction

KO validated rabbit polyclonal antibody against IDE, suitable for WB, IF/ICC, IHC-P, and ELISA. Cross-reacts with human and mouse. Ideal for studying insulin degradation, diabetes, and Alzheimer's disease.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Insulin-degrading enzyme (IDE) Rabbit pAb
Remarks/Alias	INSULYSIN
Species	Human
Gene ID	3416
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human Insulin-degrading enzyme (IDE) (NP_004960.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P14735
Protein Weight	118kDa
Shipping	Ice bag