


Ko Validated Anti-Laminin Beta 1 Rabbit Polyclonal Antibody For Wb, Elisa - P07942

Item Number: CM0024115



Introduction

Rabbit polyclonal antibody against Laminin beta 1, KO validated, suitable for WB and ELISA. Reacts with human and mouse Laminin beta 1, a key basement membrane protein involved in cell adhesion and tissue organization.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Laminin beta 1 Rabbit pAb
Remarks/Alias	CLM, LIS5, 1
Species	Human
Gene ID (Human)	3912
Gene ID	3912
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1467-1786 of human Laminin beta 1 (NP_002282.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P07942
Protein Weight	198kDa
Shipping	Ice bag