


[KO Validated] Anti-Nrbf2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q96F24

Item Number: CM0028043



Introduction

Rabbit polyclonal antibody against human NRBF2, knockout-validated, reacts with mouse and rat. Suitable for WB, IHC-P, IF/ICC, ELISA. Target involved in autophagy and transcription. Ideal for IVD core reagent research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] NRBF2 Rabbit pAb
Remarks/Alias	COPR; COPR1; COPR2; NRBF-2; F2
Species	Human
Gene ID	29982 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-287 of human NRBF2 (NP_110386.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96F24
Protein Weight	32kDa
Shipping	Ice bag