


[KO Validated] Anti-P53 Rabbit Monoclonal Antibody For Wb, Ip, If, Ihc, Chip, Elisa - P04637

Item Number: CM0011939



Introduction

High-quality rabbit monoclonal antibody against human p53 (tumor suppressor). Suitable for WB, IP, IF/ICC, IHC-P, ChIP, and ELISA. Validated KO, recognizes multiple isoforms (24-44 kDa). Ideal for cancer and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] p53 Rabbit mAb
Remarks/Alias	P53; BCC7; LFS1; BMFS5; TRP53
Species	Human
Gene ID (Human)	7157
Gene ID	7157
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-393 of human p53 (NP_000537.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ChIP, ELISA
Cross Reactivity	Human
SWISS	P04637
Protein Weight	24-44 kDa
Shipping	Ice bag