


[KO Validated] Anti-Pdk1/Pdpk1 Rabbit Pab For Wb, Ihc-P, Elisa - O15530

Item Number: CM0021039



Introduction

Rabbit polyclonal antibody targeting human PDK1/PDPK1, validated by knockout. Suitable for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. A master kinase in insulin signaling, cell proliferation, and survival.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] PDK1/PDPK1 Rabbit pAb
Remarks/Alias	PDK1; PDPK2; PDPK2P; PRO0461; K1
Species	Human
GeneID	5170 (Human)
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 409-556 of human PDK1/PDPK1 (NP_002604.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15530
Protein Weight	63kDa
Shipping	Ice bag